

MEDAPRINT FOAM[®]

Foaming disinfectant for dental impressions

MEDAPRINT FOAM is a disinfecting spray especially developed for the quick and easy cleaning and disinfection of dental impressions. Its foaming action allows a precise and easier application on the dental prints.

MEDAPRINT FOAM has been formulated to ensure a broad antimicrobial activity spectrum without affecting the dimensional accuracy of the impression materials.

It is fully compatible with alginate, silicone, polyether and hydrocolloids.

Properties

- Foam spray for easy and precise application
- Broad antimicrobial activity spectrum in 15 minutes
- Compatible with impressions made of alginate, silicone, polyether and hydrocolloids
- Effectively removes organic residues
- Bactericidal, fungicidal, tuberculocidal, mycobactericidal
- Virucidal (HBV, HIV, HCV, Herpes, Vaccinia, BVDV, influenza, Ebola, Coronavirus...)
- Does not contain phenols, aldehydes, chlorine or EDTA

Disinfecting properties

ACTIVITY SPECTRUM	STANDARD	STRAINS	CONTACT TIME
BACTERICIDAL * (Dirty conditions)	EN 13727	Pseudomonas aeruginosa Staphylococcus aureus Enterococcus hirae	5 Min.
FUNGICIDAL (Dirty conditions)	EN 13624	Candida Albicans	5 Min.
VIRUCIDAL (Dirty conditions)	DVV ⁽¹⁾ /RKI ⁽²⁾ 2014	BVDV, Vaccinia, HBV, HIV, HCV, Ebola, Herpes, Influenza H1N1, H5N1, Coronavirus	2 Min.
TUBERCULOCIDAL (Dirty conditions)	EN 14348	Mycobacterium Terrae (Surrogate M. tuberculosis)	15 Min.
MYCOBACTERICIDAL (Dirty conditions)	EN 14348	Mycobacterium Terrae Mycobacterium Avium	15 Min.

* Including all the antibiotic resistant strains as Klebsiella pneumoniae, Escherichia coli, streptococcus pneumoniae, etc.

(1) DVV: Deutsche Vereinigung zur Bekämpfung der Viruskrankheiten (German Association for the Control of Virus Diseases)

(2) RKI: Robert Koch Institute - German Federal Health Authority

Certifications

- CE mark according to the medical devices Directive (Directive 93/42/EEC)
- Medical device class IIb

For further information, please visit www.medalkan.com

MEDALKAN - Michalakopoulou 102, 115-28 Athens, Greece - Tel: (+30) 210-7484847 - Fax: (+30) 210-7772009 - www.medalkan.com - contact@medalkan.gr



Packaging

- One litre spray bottle (Ref. 20035)
- 5 litre canister (Ref. 20036)

Physical properties

- Appearance: Foaming solution
- Density: 0.99 g/cm³ at 20° C
- pH: 9.5-10.5 at 20° C
- Odour: Neutral
- Storage: 5° C - 35° C
- Stability: 3 Years
- Biodegradability: According to OCDE 301D


Instructions for use

Spray MEDAPRINT FOAM uniformly on the dental print to be disinfected to ensure a complete wetting. Let the foam act according to the desired disinfection activity (at least 2 min.). Rinse the print with running water and wipe it with a single use paper. Your dental print is clean and disinfected.

Do not mix with other disinfectants.
Store the product between + 5° C and + 35° C.

MEDAPRINT FOAM[®]

Foaming disinfectant
for dental impressions



MEDAPRINT FOAM is manufactured in the E.U.
MEDALKAN satisfies the requirements of ISO 9001:2015 for quality management system and the requirements of ISO 13485:2016 for the design and manufacture of medical devices.

Compatibility

MEDAPRINT FOAM is compatible with impressions from alginate, silicone, polyether and hydrocolloids.

MEDAPRINT FOAM is not compatible with disinfecting preparations containing aldehydes.

Composition

N- (3-aminopropyl) -N dodecylpropano-1,3-diamine,
<5% non-ionic surfactants, corrosion inhibitors, moisturising agent, excipients